

## Oximeter adapter ESAO

Ref: 452108

## **Description:**

Usable with the ESAO® 4 + interface, this adapter enables you to fully benefit from its performances: automatic air / water changeover switching, use of 3 or 4 adapters to simultaneously measure O2, CO2,

temperature, ethanol...Software or manual changeover switching of the

environment depending on the software used. Pilot lamps continuously indicate the current state of the adapter. Choice of the measurement environment Air: measurement

of the oxygen concentration in the air from 0 to 25 %.

Water: measurement of

the dissolved oxygen concentration from 0 to 20 mg /

L.Calibration

A zero potentiometer. A slope potentiometer.



## Advantages / Strengths:

Mesure précise du taux de dioxygène dans l'air et dans l'eau

## **Technical characteristics:**

Usable with the ESAO® 4 + interface, this adapter enables you to fully benefit from its performances: automatic air / water changeover switching, use of 3 or 4 adapters to simultaneously measure O2, CO2,

temperature, ethanol...Software or manual changeover switching of the environment depending on the software used.Pilot lamps

continuously indicate the current state of the adapter. Choice of the measurement environment Air: measurement

of the oxygen concentration in the air from 0 to 25 %. Water: measurement of the dissolved oxygen concentration from 0 to 20 mg / L.Calibration

A zero potentiometer. A slope potentiometer.

Mesure

 $O_2$  dans l'air (0 à 25%) et dans l'eau (0 à 20 mg/L)